Stupar Lab Meeting Multiplex CRISPR Design

Tom Kono 2017-06-19


Slides: https://z.umn.edu/crispr_slides

Publication in The Plant Cell

- Čermák et al. (2017). A Multi-purpose toolkit to enable advanced genome engineering in plants. *The Plant Cell*.
- DOI: https://doi.org/10.1105/tpc.16.00922
- See the paper for a description of the protocols, and the details of the methods
- Website: http://z.umn.edu/crisprmultiplex
 Note: Requires Javascript

Reagent Toolkit and Protocol

 A modular set of vectors that are built separately, then combined with a Golden Gate protocol



Reagent Toolkit and Protocol

Module	Number of Vectors	Description
Module A	61	Ready-to-use, with Cas9 or GFP cassettes
Module B	22	Add additional gRNAs or TALEN monomer
Module C	22	Add additional gRNAs, donor template, or expression cassettes
Transformation Backbone	31	Will be transformed into plant

Reagent Toolkit and Protocol

 Set of five protocols, with variants depending on what is being assembled (gRNAs, or TALENs)

Module B

Back to Top

Plasmid ID	Gene	Promoter	Terminator	Protocols
pMOD_B0000	None	None	None	<u>5</u>
pMOD_B2000	TALEN_2 backbone with Esp3I ccdb cassette for repeat cloning	None (begins with P2A to be fused to TALEN 1)	HSP	<u>1A,5</u>
pMOD_B2101	SapI ccdb cassette for cloning multiple gRNA protospacers with Csy4 spacers	358	35S	3A, 3S2, 5
pMOD_B2103	SapI ccdb cassette for cloning multiple gRNA protospacers with Csy4 spacers	CmYLCV	35S	3A, 3S2, 5
pMOD_B2103b	SapI ccdb cassette (promoter in the backbone, not in the assembly) for cloning multiple gRNA protospacers with Csy4 spacers	CmYLCV	35S	3A, 3S1, 3S2, 5


Vector Selection

- Use drop-down menus to choose your vector
- Menus automatically update, and it will serve a vector map as a GenBank Flat File

Vector Selection

Enter your design with the drop-down menus.

Return to main page



Primer and Map Construction

- Accepts FASTA file of target sequences
- Select target vector, promoter, restriction enzyme, and splicing system

Multi-gRNA Array Assembly - Primer Design and Map Construction

Please enter the following design parameters into the form. Note that not every vector from the Vector Selection page is available in this construction tool.

Return to the main page

Targets (FASTA):	Choose File favorite targets.fasta	
Target Vector:	pMOD_B2103	
Promoter System:	35S	\$
Restriction Enzyme:	Bsal	
Splicing System:	Csy4 ‡	
	Submit	

Primer and Map Construction

Download GenBank File

Uploaded! Your file had 6 target sequences.

You are using the 35S promoter, the BsaI resctriction enzyme, and the Csy4 splicing system. You have chosen pMOD_B2103 as your target vector.

Primer Designs

PCR Reaction 1

>o35S

TGCTCTTCGCGCATGGAGTCAAAGATTCAA

>CSY gRNA11

TGGTCTCCTGGATCTATCATCTGCCTATACGGCAGTGAAC

PCR Reaction 2

>REP gRNA11

TGGTCTCATCCAGATGTTCCGTTTTAGAGCTAGAAATAGC

>CSY gRNA13

TGGTCTCCGAAGAAGAAGAACTGCCTATACGGCAGTGAAC

PCR Reaction 3

>REP gRNA13

TGGTCTCACTTCAGACACGAGTTTTAGAGCTAGAAATAGC

>CSY_gRNA16

TGGTCTCCATATAATGCACCCTGCCTATACGGCAGTGAAC

Vector Map

```
[Download]
            pMOD B2103 favorite targets fasta
LOCUS
            MODULE B with CmYLCV: SapI ccdb cassette for cl
DEFINITION
            spacers - Csy4 .
ACCESSION
            urn
VERSION
            urn
KEYWORDS
SOURCE
  ORGANISM
                     Location/Qualifiers
FEATURES
     rep origin
                     1..857
                     /modified by="User"
                     /label="High Copy Ori"
                     878..1414
     promoter
                     /note="35S Promoter"
     misc feature
                     1422..1441
                     /note="Csy4"
     misc
                     1442..1461
                     /note="gRNA11"
     misc feature
                     1462...1537
                     /note="gRNA Repeat"
     misc feature
                     1538..1557
                     /note="Csy4"
     misc
                     1558..1577
                     /note="gRNA13"
     misc feature
```